

Anti-HDAC8 Rabbit antibody



Catalog Number: 167110

Product name

Rabbit IgG

Anti-HDAC8 Rabbit antibody

Applications

Specificity

WB, FC, IP

Human, Mouse, Rat

Dilutions

Antibody description

WB: 1/2000-1/10000

HDAC8 Rabbit polyclonal antibody

FC: 1/20

Preparation

IP: 1/20

Antigen: A synthetic peptide of human HDAC8

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage


Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Validations



Western blot detection of HDAC8 in K562 cell lysates using HDAC8 Rabbit pAb(1:1000 diluted). Predicted band size:42kDa.Observed band size:42kDa.