

Anti-DAZAP1 Rabbit antibody



Catalog Number: 167113

Product name

Anti-DAZAP1 Rabbit antibody

Specificity

Human, Rat

Antibody description

DAZAP1 Rabbit polyclonal antibody

Preparation

Antigen: Recombinant protein of human DAZAP1

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IP

Dilutions

WB: 1/2000-1/10000

IP: 1/20

Validations



Western blot detection of DAZAP1 in K562,C6,3T3,Hela cell lysates using DAZAP1 Rabbit pAb(1:1000 diluted).Predicted band size:43kDa.Observed band size:43kDa.