

# Anti-eIF4A1 Rabbit antibody



Catalog Number: 167118

## Product name

Anti-eIF4A1 Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

eIF4A1 Rabbit polyclonal antibody

## Preparation

Antigen: A synthetic peptide of human eIF4A1

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IP

## Dilutions

WB: 1/1000

IP: 1/20

## Validations



Western blot detection of eIF4A1 in K562,C6,3T3,HeLa cell lysates using eIF4A1 Rabbit pAb(1:1000 diluted).Predicted band size:46kDa.Observed band size:46kDa.