

Anti-Cleaved-Caspase-3 p17 Rabbit antibody



Catalog Number: 167126

Product name

Anti-Cleaved-Caspase-3 p17 Rabbit antibody

Ig Type

Rabbit IgG

Specificity

Human, Mouse, Rat

Applications

WB, ICC/IF, IHC-P, ELISA

Antibody description

Cleaved-Caspase-3 p17 Rabbit polyclonal antibody

Dilutions

1:500-2000

Preparation

Antigen: The antiserum was produced against synthesized peptide derived from human Caspase 3. AA range:126-175

1:50-300

1:50-300

Validations

Formulation


Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

-20°C/1 year

Clonality

Polyclonal



Western Blot analysis of 293 cells using Cleaved-Caspase-3 p17 Polyclonal Antibody