

Anti-Annexin V Rabbit pAb (HRP-Conjugate)



Catalog Number: 167129

Product name

Anti-Annexin V Rabbit pAb (HRP-Conjugate)

Specificity

Human, Monkey

Antibody description

Annexin V Rabbit pAb (HRP-Conjugate)

Preparation

Antigen: Purified recombinant human Annexin V protein fragments expressed in E.coli.

Formulation

In PBS(PH7.4) containing with 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB: 1/5000

Validations



Western blot detection of Annexin V in HeLa and COS7 cell lysates using Annexin V mouse mAb conjugated with HRP (1:5000-1:100000 diluted),incubation time:25u2103 1h.Predicted band size:30KDa.Observed band size:30KDa.