

# Anti-HDAC3 Rabbit antibody



Catalog Number: 167140

## Product name

Anti-HDAC3 Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

HDAC3 Rabbit polyclonal antibody

## Preparation

Antigen: Peptide sequence around aa.420~424(N-D-K-E-S) derived from Human HDAC3.

## Formulation

Supplied at 1.0mg/mL in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Polyclonal

## Ig Type

IgG


## Applications

WB

## Dilutions

WB:1/1000

## Validations



Western blot analysis of extracts from HeLa, Molt-4 and 3T3 cells using HDAC3 rabbit pAb (dilution 1:1000). Predicted band size: 49kDa. Observed band size: 49kDa.