

Anti-HDAC3 Rabbit antibody



Catalog Number: 167140

Product name

Anti-HDAC3 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

HDAC3 Rabbit polyclonal antibody

Preparation

Antigen: Peptide sequence around aa.420~424(N-D-K-E-S) derived from Human HDAC3.

Formulation

Supplied at 1.0mg/mL in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB:1/1000

Validations



Western blot analysis of extracts from HeLa, Molt-4 and 3T3 cells using HDAC3 rabbit pAb (dilution 1:1000). Predicted band size: 49 kDa. Observed band size: 49 kDa.