

# Anti-P4HB Rabbit antibody



Catalog Number: 167161

## Product name

Anti-P4HB Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

P4HB Rabbit polyclonal antibody

## Preparation

Antigen: A synthetic peptide of human PDIA2

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, FC

## Dilutions

WB: 1/2000-1/10000

FC: 1/20

## Validations



Western blot detection of PDIA2 in K562,C6 cell lysates using PDIA2 Rabbit pAb(1:1000 diluted).Predicted band size:58kDa.Observed band size:58kDa.