

Anti-Bag3 Rabbit antibody

Catalog Number: 167162

Product name

Anti-Bag3 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

Bag3 Rabbit polyclonal antibody

Preparation

Antigen: Recombinant protein of human Bag3

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, FC

Dilutions

WB: 1/1000

FC: 1/20

Validations



Western blot detection of Bag3 in K562,Rat Brain cell lysates using Bag3 Rabbit pAb(1:1000 diluted).Predicted band size:62kDa.Observed band size:80kDa.