

Anti-PDGFR beta Rabbit antibody



Catalog Number: 167163

Product name

Anti-PDGFR beta Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

PDGFR beta Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human PDGFR beta

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M Nacl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, ICC/IF, FC, IP

Dilutions


WB: 1/1000

ICC/IF: 1/20

FC: 1/20

IP: 1/20

Validations



Western blot detection of PDGFR beta in 3T3 cell lysates using PDGFR beta Rabbit pAb(1:1000 diluted).Predicted band size:124KDa.Observed band size:190KDa.