

# Anti-Tec Rabbit antibody

Catalog Number: 167168

## Product name

Anti-Tec Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Tec Rabbit polyclonal antibody

## Preparation

Antigen: A synthetic peptide of human Tec

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, ICC/IF, FC


## Dilutions

WB: 1/2000


ICC/IF: 1/20-1/100

FC: 1/20

## Validations



Immunocytochemistry of Tec(green) in HeLa cells using Tec Rabbit pAb at dilution 1/50, and



Western blot detection of Tec in K562 cell lysates using Tec Rabbit pAb(1:1000 diluted). Predicted band size:74kDa.Observed band size:58-75kDa.