

Anti-Dnmt1 Rabbit antibody

Catalog Number: 167177

Product name

Anti-Dnmt1 Rabbit antibody

Specificity

Human

Antibody description

Dnmt1 Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human Dnmt1

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, ICC/IF, FC, IP

Dilutions


WB: 1/1000

ICC/IF: 1/20


FC: 1/20-1/100

IP: 1/20

Validations



Immunocytochemistry of Dnmt1(green) in HeLa cells using Dnmt1 Rabbit pAb at dilution 1/50, and



DAPI(blue)

Western blot detection of Dnmt1 in Jurkat cell lysates using Dnmt1 Rabbit pAb(1:1000 diluted).Predicted band size:183KDa.Observed band size:183KDa.