

Anti-K-48 Linkage Specific Ubiquitin Rabbit antibody



Catalog Number: 167195

Product name

Anti-K-48 Linkage Specific Ubiquitin Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

K-48 Linkage Specific Ubiquitin Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human K-48 Linkage Specific Ubiquitin

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, FC

Dilutions


WB: 1/2000

IHC: 1/50

ICC/IF: 1/100

FC: 1/20

Validations



Immunohistochemistry of K-48 Linkage Specific Ubiquitin in paraffin-embedded Human colon cancer tissue using K-48 Linkage Specific

Ubiquitin Rabbit pAb at dilution 1/50

Western blot detection of K-48 Linkage Specific Ubiquitin in 293 cell lysates using K-48 Linkage Specific Ubiquitin Rabbit pAb(1:1000 diluted).

