

Anti-Thymidylate Synthase Rabbit antibody



Catalog Number: 167226

Product name

Anti-Thymidylate Synthase Rabbit antibody

WB: 1/1000-1/5000

IHC: 1/20-1/50

Specificity

Human, Mouse, Rat

ICC/IF: 1/20-1/50

FC: 1/20

IP: 1/20-1/50

Antibody description

Thymidylate Synthase Rabbit polyclonal antibody

Validations

Preparation

Antigen: A synthetic peptide of human Thymidylate Synthase

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal


Ig Type

Rabbit IgG


Applications

WB, IHC-P, ICC/IF, FC, IP

Dilutions



Immunohistochemistry of Thymidylate Synthase in paraffin-embedded Human tonsil using Thymidylate Synthase Rabbit pAb at dilution 1/20



Western blot detection of Thymidylate Synthase in K562 cell lysates using Thymidylate Synthase Rabbit pAb(1:1000 diluted).Predicted band size:36kDa.Observed band size:36kDa.