

Anti-VCP Rabbit antibody



Catalog Number: 167237

Product name

Anti-VCP Rabbit antibody

IHC: 1/20-1/100

ICC/IF: 1/20-1/50

Specificity

Human, Mouse, Rat

FC: 1/20

IP: 1/20-1/50


Antibody description

VCP Rabbit polyclonal antibody

Validations

Preparation

Antigen: A synthetic peptide of human VCP



Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Western blot detection of VCP in MCF7,C6,3T3 cell lysates using VCP Rabbit pAb(1:1000 diluted).Predicted band size:89KDa.Observed band size:89KDa.

Clonality


Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, FC, IP



Dilutions

WB: 1/2000-1/10000

Immunohistochemistry of VCP in paraffin-embedded Human tonsil using VCP Rabbit pAb at dilution 1/20