

# Anti-TMEM192 Rabbit antibody



Catalog Number: 167243

## Product name

Anti-TMEM192 Rabbit antibody

## Specificity

Human

## Antibody description

TMEM192 Rabbit polyclonal antibody

## Preparation

Antigen: Recombinant protein of human TMEM192

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, FC

## Dilutions

WB: 1/1000-1/2000

FC: 1/20

## Validations



Western blot detection of TMEM192 in K562 cell lysates using TMEM192 Rabbit pAb(1:1000 diluted). Predicted band size:31kDa.Observed band size:31kDa.