

# Anti-PRKACG Rabbit antibody



Catalog Number: 167245

## Product name

Anti-PRKACG Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

PRKACG Rabbit polyclonal antibody

## Preparation

Antigen: Recombinant protein of human PRKACG

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, FC, IP


## Dilutions

WB: 1/5000-1/10000

FC: 1/20-1/100

IP: 1/20-1/50

## Validations



Western blot detection of PRKACG in MCF7,C6,HeLa cell lysates using PRKACG Rabbit pAb(1:1000 diluted).Predicted band size:40KDa.Observed band size:40KDa.