

# Anti-Coilin Rabbit antibody



Catalog Number: 167250

## Product name

Anti-Coilin Rabbit antibody

## Specificity

Human

## Antibody description

Coilin Rabbit polyclonal antibody

## Preparation

Antigen: A synthetic peptide of human Coilin

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, FC

## Dilutions

WB: 1/1000

FC: 1/50

## Validations



Western blot detection of Coilin in K562 cell lysates using Coilin Rabbit pAb(1:1000 diluted).Predicted band size:63kDa.Observed band size:80kDa.