

Anti-Cdk4 Rabbit antibody



Catalog Number: 167253

General Information

Product name:

Anti-Cdk4 Rabbit antibody

Immunogen:

Applications:

WB, ICC/IF, FC

Specificity:

Human

Antibody description:

Cdk4 Rabbit polyclonal antibody

Properties

Preparation:

Antigen: A synthetic peptide of human Cdk4

Formulation:

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage:

Store at -20°C. Stable for 12 months from date of receipt.

Clonality:

Polyclonal

Ig Type:Rabbit IgG


Applications

WB: 1/1000-1/5000

ICC/IF: 1/100

FC: 1/20

Validations



Western blot detection of Cdk4 in HeLa cell lysates using Cdk4 Rabbit pAb(1:1000 diluted).Predicted band size:34kDa.Observed band size:34kDa.