

Anti-Cdk4 Rabbit antibody



Catalog Number: 167259

General Information

Product name:

Anti-Cdk4 Rabbit antibody

Immunogen:

Applications:

WB, ICC/IF, FC, IP

Specificity:

Human, Mouse, Rat

Antibody description:

Cdk4 Rabbit polyclonal antibody

Properties

Preparation:

Antigen: A synthetic peptide of human Cdk4

Formulation:

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage:

Store at -20°C. Stable for 12 months from date of

receipt.

Clonality:

Polyclonal

Ig Type:Rabbit IgG

Applications


WB: 1/1000

ICC/IF: 1/50

FC: 1/20

IP: 1/20

Validations



Western blot detection of Cdk4 in K562, HeLa cell lysates using Cdk4 Rabbit pAb(1:1000 diluted). Predicted band size: 34kDa. Observed band size: 34kDa.