

Anti-ATG7 Rabbit antibody



Catalog Number: 167266

Product name

Anti-ATG7 Rabbit antibody

Specificity

Human

Antibody description

ATG7 Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human ATG7

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, ICC/IF, FC, IP

Dilutions


WB: 1/2000-1/10000

ICC/IF: 1/20-1/100

FC: 1/20

IP: 1/20

Validations



Western blot detection of ATG7 in K562 cell lysates using ATG7 Rabbit pAb(1:1000 diluted). Predicted band size:78kDa. Observed band size:78kDa.