

Anti-Ubiquitin Rabbit antibody



Catalog Number: 167270

Product name

Anti-Ubiquitin Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

Ubiquitin Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human Ubiquitin

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, ICC/IF, FC


Dilutions

WB: 1/1000

ICC/IF: 1/50

FC: 1/20-1/100

Validations



Western blot detection of Ubiquitin in 293 cell lysates using Ubiquitin Rabbit pAb(1:1000 diluted).Predicted band size:8KDa.Observed band size:8KDa.