

# Anti-MLH1 Rabbit antibody



Catalog Number: 167278

## General Information

### Product name:

Anti-MLH1 Rabbit antibody

### Immunogen:

### Applications:

WB, IHC-P, ICC/IF, FC, IP

### Specificity:

Human

### Antibody description:

MLH1 Rabbit polyclonal antibody

## Properties

### Preparation:

Antigen: A synthetic peptide of human MLH1

### Formulation:

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

### Storage:

Store at -20°C. Stable for 12 months from date of receipt.

### Clonality:

Polyclonal antibody

**Ig Type:**Rabbit IgG

## Applications

WB: 1/1000


IHC: 1/20

ICC/IF: 1/20

FC: 1/100

IP: 1/20

## Validations



Western blot detection of MLH1 in K562,C6,3T3,HeLa cell lysates using MLH1 Rabbit pAb(1:1000 diluted).Predicted band size:85kDa.Observed band size:85kDa.