

Anti-TRK fused gene Rabbit antibody



Catalog Number: 167281

Product name

Anti-TRK fused gene Rabbit antibody

WB: 1/1000-1/5000

IHC: 1/50

Specificity

Human, Mouse, Rat

ICC/IF: 1/100

FC: 1/20

Antibody description

TRK fused gene Rabbit polyclonal antibody

Validations

Preparation

Antigen: A synthetic peptide of human TRK fused gene

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal


Ig Type

Rabbit IgG


Applications

WB, IHC-P, ICC/IF, FC

Dilutions



Immunohistochemistry of TRK fused gene in paraffin-embedded Human tonsil using TRK fused



gene Rabbit pAb at dilution 1/50

Western blot detection of TRK fused gene in C6,3T3 cell lysates using TRK fused gene Rabbit pAb(1:1000 diluted).Predicted band size:44kDa.Observed band size:44kDa.