

# Anti-Peroxiredoxin 3 Rabbit antibody



Catalog Number: 167301

## Product name

Anti-Peroxiredoxin 3 Rabbit antibody

## Specificity

Human

## Antibody description

Peroxiredoxin 3 Rabbit polyclonal antibody

## Preparation

Antigen: A synthesized peptide derived from human Peroxiredoxin 3

## Formulation

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P, ICC/IF, IP, FC

## Dilutions

WB: 1/1000-1/5000


IHC: 1/50-1/200

ICC/IF: 1/50-1/200

IP: 1/80

FC: 1/50

## Validations



Western blot detection of Peroxiredoxin 3 in K562,Raji,MCF7 cell lysates using Peroxiredoxin 3 Rabbit pAb(1:1000 diluted).Predicted band size:26kDa.Observed band size:26kDa.