

# Anti-Hsp70 Rabbit antibody



Catalog Number: 167324

## Product name

Anti-Hsp70 Rabbit antibody

## Specificity

Human

## Antibody description

Hsp70 Rabbit polyclonal antibody

## Preparation

Antigen: Recombinant protein of human Hsp70

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P, ICC/IF, FC, IP

## Dilutions

WB: 1/1000


IHC: 1/50

ICC/IF: 1/20

FC: 1/20

IP: 1/20

## Validations



Western blot detection of Hsp70 in HeLa, HCT116, A549 cell lysates using Hsp70 Rabbit pAb(1:1000 diluted). Predicted band size: 70KDa. Observed band size: 70KDa.