

Anti-CNOT7 Rabbit antibody



Catalog Number: 167325

Product name

Anti-CNOT7 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

CNOT7 Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of mouse CNOT7

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, FC, IP

Dilutions


WB: 1/1000

IHC: 1/200

FC: 1/100

IP: 1/20

Validations



Western blot detection of CNOT7 in Molt-4,Raji,COS7,Hela cell lysates using CNOT7 Rabbit pAb(1:1000 diluted).Predicted band size:33KDa.Observed band size:33KDa.