

Anti-eIF4EBP1 Rabbit antibody



Catalog Number: 167326

Product name

Anti-eIF4EBP1 Rabbit antibody

Specificity

Human

Antibody description

eIF4EBP1 Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human eIF4EBP1

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, FC, IP

Dilutions


WB: 1/1000

IHC: 1/20

FC: 1/20

IP: 1/20

Validations



Western blot detection of 4E-BP1 in K562,293,HeLa cell lysates using 4E-BP1 Rabbit pAb(1:1000 diluted).Predicted band size:12KDa.Observed band size:15-20KDa.