

Anti-ELMO1 Rabbit antibody



Catalog Number: 167345

Product name

Anti-ELMO1 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

ELMO1 Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human ELMO1

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, ICC/IF, FC, IP

Dilutions


WB: 1/5000

ICC/IF: 1/50

FC: 1/20

IP: 1/20

Validations



Western blot detection of ELMO1 in Rat Brain,C6,3T3 cell lysates using ELMO1 Rabbit pAb(1:1000 diluted).Predicted band size:84kDa.Observed band size:84kDa.