

Anti-CXCR5 Rabbit antibody



Catalog Number: 167346

Product name

Anti-CXCR5 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

CXCR5 Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human CXCR5

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, ICC/IF, FC


Dilutions

WB: 1/1000-1/5000

ICC/IF: 1/50

FC: 1/20

Validations



Western blot detection of CXCR5 in K562,C6,3T3,Hela cell lysates using CXCR5 Rabbit pAb(1:1000 diluted).Predicted band size:42kDa.Observed band size:42kDa.