

Anti-AGPS Rabbit antibody



Catalog Number: 167352

Product name

Anti-AGPS Rabbit antibody

WB: 1/1000-1/5000

IHC: 1/20

Specificity

Human, Mouse, Rat

FC: 1/20

IP: 1/20

Antibody description

AGPS Rabbit polyclonal antibody

Validations

Preparation

Antigen: Recombinant protein of human Alkyl-DHAP synthase

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal


Ig Type

Rabbit IgG


Applications

WB, IHC-P, FC, IP

Dilutions



Western blot detection of Alkyl-DHAP synthase in K562 cell lysates using Alkyl-DHAP synthase Rabbit pAb(1:1000 diluted). Predicted band size: 73kDa. Observed band size: 73kDa.



Immunohistochemistry of Alkyl-DHAP synthase in paraffin-embedded Human colon cancer tissue using Alkyl-DHAP synthase Rabbit pAb at dilution 1/20