

Anti-Cytokeratin 7 Rabbit antibody



Catalog Number: 167385

Product name

Anti-Cytokeratin 7 Rabbit antibody

ICC/IF: 1/50-1/200

IP: 1/20

Specificity

Human

Validations

Antibody description

Cytokeratin 7 Rabbit polyclonal antibody

Preparation

Antigen: A synthesized peptide derived from human Cytokeratin 7

Formulation

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB, IHC-P, ICC/IF, IP, FC


Dilutions

WB: 1/5000-1/10000

IHC: 1/50-1/200




Immunohistochemical analysis of paraffin-embedded human breast carcinoma, using



Cytokeratin 7 Antibody .

Western blot analysis of Cytokeratin 7 expression



in T47 D cell lysate.

Immunofluorescent analysis of HeLa cells, using Cytokeratin 7 Antibody.