

# Anti-Cytokeratin 7 Rabbit antibody



Catalog Number: 167385

## Product name

Anti-Cytokeratin 7 Rabbit antibody


ICC/IF: 1/50-1/200

IP: 1/20

## Validations

## Specificity

Human



Immunohistochemical analysis of paraffin-embedded human breast carcinoma, using

## Antibody description

Cytokeratin 7 Rabbit polyclonal antibody

## Preparation

Antigen: A synthesized peptide derived from human Cytokeratin 7

## Formulation

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB, IHC-P, ICC/IF, IP, FC

## Dilutions

WB: 1/5000-1/10000


IHC: 1/50-1/200



Cytokeratin 7 Antibody .

Western blot analysis of Cytokeratin 7 expression

in T47 D cell lysate.



Immunofluorescent analysis of HeLa cells, using Cytokeratin 7 Antibody.