

Anti-Fascin Rabbit antibody



Catalog Number: 167409

Product name

Anti-Fascin Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

Fascin Rabbit polyclonal antibody

Preparation

Antigen: Recombinant protein of human Fascin

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, FC, IP

Dilutions

WB: 1/2000-1/10000


IHC: 1/50-1/200

ICC/IF: 1/50

FC: 1/20-1/100

IP: 1/20-1/50

Validations



Western blot detection of Fascin in HeLa, A549, HepG2, COS7 cell lysates using Fascin Rabbit pAb(1:1000 diluted). Predicted band size: 54KDa. Observed band size: 54KDa.