

# Anti-NFkB p100/p52 Rabbit antibody



Catalog Number: 167418

## Product name

Anti-NFkB p100/p52 Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

NFkB p100/p52 Rabbit polyclonal antibody

## Preparation

Antigen: A synthesized peptide derived from human NFkB p100 / p52

## Formulation

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, ICC/IF, IP


## Dilutions

WB: 1/1000-1/5000

ICC/IF: 1/50-1/200

IP: 1/50

## Validations



Western blot analysis of NFkB p100 / p52 expression in (1) Jurkat cell lysate; (2) Mouse heart lysate.