

Anti-LGALS14 Rabbit antibody



Catalog Number: 167423

Product name

Anti-LGALS14 Rabbit antibody

Specificity

Human, Mouse

Antibody description

LGALS14 Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human LGALS14

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, ICC/IF

Dilutions

WB: 1/1000

ICC/IF: 1/50

Validations



Western blot detection of LGALS14 in C6 cell lysates using LGALS14 Rabbit pAb(1:1000 diluted). Predicted band size:16kDa. Observed band size:16kDa.