

Anti-NDUFAF4 Rabbit antibody



Catalog Number: 167449

Product name

Anti-NDUFAF4 Rabbit antibody

Specificity

Human

Antibody description

NDUFAF4 Rabbit polyclonal antibody

Preparation

Antigen: Recombinant protein of human NDUFAF4

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, ICC/IF, IP


Dilutions

WB: 1/1000

ICC/IF: 1/20

IP: 1/20

Validations



Western blot detection of NDUFAF4 in K562, HeLa cell lysates using NDUFAF4 Rabbit pAb(1:1000 diluted). Predicted band size: 20kDa. Observed band size: 20kDa.