

Anti-RAB22A Rabbit antibody



Catalog Number: 167477

Product name

Anti-RAB22A Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

RAB22A Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human RAB22A

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, ICC/IF, IP


Dilutions

WB: 1/1000

ICC/IF: 1/50

IP: 1/20

Validations



Western blot detection of RAB22A in Rat Brain,C6,3T3,Hela cell lysates using RAB22A Rabbit pAb(1:1000 diluted).Predicted band size:22kDa.Observed band size:22kDa.