

# Anti-NTPCR Rabbit antibody



Catalog Number: 167479

## Product name

Anti-NTPCR Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

NTPCR Rabbit polyclonal antibody

## Preparation

Antigen: Recombinant protein of human Nucleoside triphosphate phosphohydrolase

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, ICC/IF, IP


## Dilutions

WB: 1/1000

ICC/IF: 1/20

IP: 1/20

## Validations



Western blot detection of Nucleoside triphosphate phosphohydrolase in K562, C6, 3T3, HeLa cell lysates using Nucleoside triphosphate phosphohydrolase Rabbit pAb(1:1000 diluted). Predicted band size: 21 kDa. Observed band size: 21 kDa.