

# Anti-Histone H1.0 Rabbit antibody



Catalog Number: 167482

## Product name

Anti-Histone H1.0 Rabbit antibody

## Specificity

Human

## Antibody description

Histone H1.0 Rabbit polyclonal antibody

## Preparation

Antigen: A synthetic peptide of human Histone H1.0

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, ICC/IF

## Dilutions

WB: 1/1000-1/5000

ICC/IF: 1/20-1/50

## Validations



Western blot detection of Histone H1.0 in HeLa cell lysates using Histone H1.0 Rabbit pAb(1:1000 diluted). Predicted band size:21kDa.Observed band size:30kDa.