

Anti-BNIP3L Rabbit antibody



Catalog Number: 167489

Product name

Anti-BNIP3L Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

BNIP3L Rabbit polyclonal antibody

Preparation

Antigen: A synthesized peptide derived from human BNIP3L

Formulation

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC-P

Dilutions

WB: 1/500-1/2000

IHC: 1/50-1/200

Validations



Western blot analysis of BNIP3L expression in K562 cell lysate.