

# Anti-Iba1 Rabbit antibody

Catalog Number: 167492

## Product name

Anti-Iba1 Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Iba1 Rabbit polyclonal antibody

## Preparation

Antigen: A synthetic peptide of human Iba1

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IHC-P

## Dilutions


WB: 1/1000

IHC: 1/20-1/500

## Validations



Western blot detection of Iba1 in C6,3T3,Hela cell lysates using Iba1 Rabbit pAb(1:1000 diluted).Predicted band size:17kDa.Observed band size:17kDa.



Immunohistochemistry of Iba1 in paraffin-embedded Human colon cancer tissue using Iba1 Rabbit pAb at dilution 1/20