

# Anti-Tetranectin Rabbit antibody



Catalog Number: 167503

## Product name

Anti-Tetranectin Rabbit antibody

WB: 1/1000-1/5000

IHC: 1/50-1/200

## Specificity

Human, Mouse, Rat

IP: 1/20-1/50

## Antibody description

Tetranectin Rabbit polyclonal antibody

## Preparation

Antigen: A synthetic peptide of human Tetranectin

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type


Rabbit IgG

## Applications


WB, IHC-P, IP

## Dilutions

## Validations



Western blot detection of Tetranectin in MCF7 cell lysates using Tetranectin Rabbit pAb(1:1000 diluted). Predicted band size:23KDa. Observed band size:23KDa.



Immunohistochemistry of Tetranectin in paraffin-embedded Human lung cancer tissue using Tetranectin Rabbit pAb at dilution 1/50