

Anti-MOG Rabbit antibody



Catalog Number: 167504

Product name

Anti-MOG Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

MOG Rabbit polyclonal antibody

Preparation

Antigen: Recombinant protein of human Myelin oligodendrocyte glycoprotein

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, IP


Dilutions

WB: 1/1000

IHC: 1/100

IP: 1/20

Validations



Western blot detection of Myelin oligodendrocyte glycoprotein in C6 cell lysates using Myelin oligodendrocyte glycoprotein Rabbit pAb(1:1000 diluted). Predicted band size:28KDa. Observed band size:28KDa.