

# Anti-Artemin Rabbit antibody



Catalog Number: 167506

## Product name

Anti-Artemin Rabbit antibody

WB: 1/1000

IHC: 1/20

## Specificity

Human, Mouse, Rat

IP: 1/20

## Antibody description

Artemin Rabbit polyclonal antibody

## Preparation

Antigen: A synthetic peptide of human Artemin

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type


Rabbit IgG

## Applications


WB, IHC-P, IP

## Dilutions

## Validations



Immunohistochemistry of Artemin in paraffin-embedded Human lung cancer tissue using Artemin Rabbit pAb at dilution 1/20



Western blot detection of Artemin in K562,Rat Brain,C6,3T3,Hela cell lysates using Artemin Rabbit pAb(1:1000 diluted).Predicted band size:23kDa.Observed band size:23kDa.