

Anti-Phospho-Glycogen synthase (Ser641) Rabbit antibody



Catalog Number: 167509

Product name

Anti-Phospho-Glycogen synthase (Ser641) Rabbit antibody

Specificity

Human, Mouse

Antibody description

Phospho-Glycogen synthase (Ser641) Rabbit polyclonal antibody

Preparation

Antigen: A synthetic phosphopeptide corresponding to residues surrounding Ser641 of human Glycogen synthase 1/GYS1

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, IP


Dilutions

WB: 1/1000

IHC: 1/50

IP: 1/20

Validations



Western blot detection of Glycogen synthase 1/GYS1 (Phospho-Ser641) in C6,3T3 cell lysates using Glycogen synthase 1/GYS1 (Phospho-Ser641) Rabbit pAb(1:1000 diluted). Predicted band size:84kDa.Observed band size:84kDa.