

# Anti-Histone H2B Rabbit antibody



Catalog Number: 167513

## Product name

Anti-Histone H2B Rabbit antibody

## Specificity

Human, Mouse

## Antibody description

Histone H2B Rabbit polyclonal antibody

## Preparation

Antigen: A synthetic peptide of human Histone H2B

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P, IP


## Dilutions

WB: 1/2000

IHC: 1/200

IP: 1/20

## Validations



Western blot detection of Histone H2B in K562,C6 cell lysates using Histone H2B Rabbit pAb(1:1000 diluted).Predicted band size:14kDa.Observed band size:14kDa.