

# Anti-A1BG Rabbit antibody



Catalog Number: 167516

## Product name

Anti-A1BG Rabbit antibody

## Specificity

Human

## Antibody description

A1BG Rabbit polyclonal antibody

## Preparation

Antigen: Synthesized peptide derived from A1BG .  
at AA range: 100-180

## Formulation

Liquid in PBS containing 50% glycerol, 0.5% BSA  
and 0.02% sodium azide.

## Storage

Store at 4°C short term. Store at -20°C long term.  
Avoid freeze / thaw cycle.

## Clonality

Polyclonal

## Ig Type

Unknown

## Applications

WB, IHC-P, ELISA


## Dilutions

WB: 1/500 - 1/2000


IHC-p: 1/100-300

ELISA: 1/40000


## Validations



Western Blot analysis of extracts from K562 cells, using A1BG Polyclonal Antibody. Secondary antibody was diluted at 1:20000



Immunohistochemical analysis of paraffin-embedded human-uterus, antibody was diluted at




1:100

Immunohistochemical analysis of paraffin-embedded human-ovary, antibody was diluted at

# Anti-A1BG Rabbit antibody


Catalog Number: 167516

---



1:100

Immunohistochemical analysis of paraffin-embedded human-uterus, antibody was diluted at



1:100

Immunohistochemical analysis of paraffin-embedded human-ovary, antibody was diluted at 1:100