

# Anti-MURF1 Rabbit antibody



Catalog Number: 167517

## Product name

Anti-MURF1 Rabbit antibody

;IHC: 1/50 - 1/200

;IF: 1/50 - 1/200

## Specificity

Human, Mouse, Rat

## Antibody description

MURF1 Rabbit polyclonal antibody

## Preparation

Antigen: A synthetic peptide of human TRIM63

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C . Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB, IHC-P, ICC/IF


## Dilutions

WB: 1/500 - 1/2000

## Validations



Immunofluorescence - TRIM63 Polyclonal Antibody. Immunofluorescence analysis of A549 cells using TRIM63 antibody . Blue: DAPI for



nuclear staining.

Western blot - TRIM63 Polyclonal Antibody. Western blot analysis of extracts of mouse heart, using TRIM63 antibody at 1:400 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit . Exposure time: 90s.