

# Anti-Phospho-Rb (Ser780) Rabbit antibody



Catalog Number: 167528

## Product name

Anti-Phospho-Rb (Ser780) Rabbit antibody

## Ig Type

Rabbit IgG

## Specificity

Human

## Applications

WB, IHC-P

## Antibody description

Phospho-Rb (Ser780) Rabbit polyclonal antibody

## Dilutions

WB: 1/500-1/2000

## Preparation

Antigen: A synthesized peptide derived from human Retinoblastoma

IHC: 1/50-1/200

## Formulation

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.


## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Polyclonal

## Validations



Western blot analysis of Retinoblastoma phosphorylation expression in Jurkat cell lysate treated with Alkaline Phosphatase.